

SETBP1 antibody



Catalog Number: 115220

Product name

SETBP1 antibody

Immunogen

[Human SETBP1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SETBP1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 145-150 kDa

Preparation

This antibody was obtained by immunization of SETBP1 recombinant protein (Accession Number: NM_001130110). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115220(SETBP1 Antibody) at dilution of 1:300