

# Flag tag antibody



Catalog Number: 107220

## Product name

Flag tag antibody

## Specificity

Human; other species not tested.

## Antibody description

Flag tag Mouse Monoclonal antibody. Positive IF detected in Transfected HEK-293 cells. Positive WB detected in Recombinant protein, Transfected HEK-293 cells. Positive IP detected in Recombinant protein.

## Preparation

This antibody was obtained by immunization of Peptide. Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications

ELISA, WB, IP, IF


## Dilutions

Recommended Dilution:


WB: 1:5000-1:50000

IP: N/AIF: 1:20-1:200


## Validations



IP Result of anti-Flag tag (IP:Catalog No:117312, 3ug; Detection:Catalog No:117313 1:2000) with Recombinant protein protein lysate 1280ug.



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117312(Flag tag Antibody) at dilution of 1:1000




Recombinant protein were Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:117312 (Flag tag Antibody) at various dilutions.

# Flag tag antibody



Catalog Number: 107220

---



Immunofluorescent analysis of Transfected HEK-293 cells using Catalog No:117312(Flag tag Antibody) at dilution of 1:50 and Alexa Fluor 594-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)