

GNA13 Polyclonal Antibody



Catalog Number: 164620

Product name

GNA13 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to GNA13

Preparation

Antigen: Recombinant protein of human GNA13

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:100 - 1:200


IF 1:50 - 1:200

Validations



Western blot - GNA13 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GNA13 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 90s.


Western blot - GNA13 Polyclonal Antibody

Western blot analysis of extracts from normal (control) and GNA13 knockout HeLa cells, using GNA13 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time:

GNA13 Polyclonal Antibody



Catalog Number: 164620




10s.

Immunohistochemistry - GNA13 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using GNA13 antibody (40x


human esophagus using GNA13 antibody (40x



lens).

Immunohistochemistry - GNA13 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using GNA13 antibody (40x lens).



lens).

Immunohistochemistry - GNA13 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded

Immunohistochemistry - GNA13 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using GNA13 antibody (40x lens).