

RLBP1L2 antibody



Catalog Number: 114720

Product name

RLBP1L2 antibody

Immunogen

[Human RLBP1L2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

RLBP1L2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, Jurkat cells, L02 cells. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 38kd

Preparation

This antibody was obtained by immunization of RLBP1L2 recombinant protein (Accession Number: NM_001010852). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114720(RLBP1L2 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using RLBP1L2 antibody Catalog No:114720 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).