

HLA-DRB3 Polyclonal Antibody



Catalog Number: 163030

Product name

HLA-DRB3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to HLA-DRB3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 30-227 of human HLA-DRB3 (NP_072049.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - HLA-DRB3 Polyclonal Antibody

Western blot analysis of extracts of Raji cells, using HLA-DRB3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunohistochemistry - HLA-DRB3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded tonsil using HLA-DRB3 antibody at dilution of



1:100 (40x lens).


HLA-DRB3 Polyclonal Antibody



Catalog Number: 163030

Immunohistochemistry - HLA-DRB3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using HLA-DRB3 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - HLA-DRB3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human appendix using HLA-DRB3 antibody at dilution of 1:100 (40x lens).