

# TSGA10 antibody



Catalog Number: 116430

## Product name

TSGA10 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

TSGA10 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse brain tissue. Observed molecular weight by Western-blot: 81kd

## Preparation

This antibody was obtained by immunization of TSGA10 recombinant protein (Accession Number: XM\_017005039). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116430(TSGA10 antibody) at dilution of 1:300