

Anti-BRAT1 Rabbit antibody



Catalog Number: 166830

Product name

Anti-BRAT1 Rabbit antibody

Specificity

Human

Antibody description

BRAT1 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human BRAT1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/2000-1/10000

ICC/IF: 1/50


FC: 1/20

IP: 1/20

Validations



Immunocytochemistry of BRAT1(green) in Hela cells using BRAT1 Rabbit pAb at dilution 1/50, and



Western blot detection of BRAT1 in K562 cell lysates using BRAT1 Rabbit pAb(1:1000 diluted). Predicted band size:88kDa. Observed band size:88kDa.