

DDX4,VASA antibody



Catalog Number: 109830

Product name

DDX4,VASA antibody

Immunogen

[Human DDX4,VASA Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

DDX4,VASA Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 79 kDa

Preparation

This antibody was obtained by immunization of DDX4,VASA recombinant protein (Accession Number: NM_001142549). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109830(DDX4,VASA antibody) at dilution of



1:200

IP Result of anti-DDX4,VASA (IP:Catalog No:109830, 4ug; Detection:Catalog No:109830 1:300) with mouse testis tissue lysate 4000ug.