

Dynamitin antibody



Catalog Number: 110130

Product name

Dynamitin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Dynamitin Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, human brain tissue, human heart tissue, human testis tissue. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of Dynamitin recombinant protein (Accession Number: AY189155). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110130(dynactin-2 antibody) at dilution of 1:500