

# Recombinant PSB5 protein (His tag)



Catalog Number: 544530

---

## Gene Name Synonym

Macropain epsilon chain; Multicatalytic endopeptidase complex epsilon chain; Proteasome chain 6; Proteasome epsilon chain; Proteasome subunit MB1; Proteasome subunit X; Proteasome subunit beta-5; beta-5

## Protein Construction

DNA encoding 60-263 aa (P28074) were fused with 6His tag.

## Organism

Human

## Expression Host

E. coli

## Activity

Not tested.

## Purity

>86%, by SDS-PAGE with Coomassie Brilliant Blue staining.

## Endotoxin

Please contact the lab for more information.

## Stability

Samples are stable for up to twelve months from date of receipt at -70°C

## Molecular Mass

Predictes a molecular mass of 24.11 kDa. In SDS-PAGE under reducing conditions, it migrates as an approximately 24 kDa band.

## Formulation

Lyophilized from sterile PBS, pH 7.4

## Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

## Reconstitution

Reconstitute at 0.25 µg/µl in sterile water for short-term storage. Reconstitution with 50% glycerol solution is recommended for longer term storage (see Stability and Storage for more details). If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used). Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.