

# C1orf104 antibody



Catalog Number: 108630

## Product name

C1orf104 antibody

## Immunogen

[Human C1orf104 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

C1orf104 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 30 kDa

## Preparation

This antibody was obtained by immunization of C1orf104 recombinant protein (Accession Number: BC080538). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108630(C1orf104 Antibody) at dilution of 1:1000