

HACE1 antibody



Catalog Number: 111330

Product name

HACE1 antibody

Immunogen

[Human HACE1 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

HACE1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 78 kDa

Preparation

This antibody was obtained by immunization of HACE1 recombinant protein (Accession Number: NM_020771). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111330(HACE1 Antibody) at dilution of 1:600



Immunofluorescent analysis of HeLa cells using Catalog No:111330(HACE1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)