

Transferrin antibody



Catalog Number: 116330

Product name

Transferrin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Transferrin Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, rat kidney tissue. Positive IHC detected in human hepat°Cirrhosis tissue, human placenta tissue, human testis tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of Transferrin recombinant protein (Accession Number: XM_017007089). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:117341(TF antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:117341(TF antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:117341(TF antibody) at dilution of 1:50 (under 10x lens)

Transferrin antibody



Catalog Number: 116330

40x lens)