

# Anti-CHD4 (5F1) Mouse antibody

Catalog Number: 168330

## Product name

Anti-CHD4 (5F1) Mouse antibody

## Specificity

Human

## Antibody description

CHD4 (5F1) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human CHD4 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in PBS(pH 7.4)containing with 0.02% sodium azide 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

WB, IP


## Dilutions

WB: 1/1000

## Validations



Western blot detection of CHD4 in K562, HeLaNE, HL-60, MOLT-4, Jurkat and 293 cell lysates using CHD4 mouse mAb (1:1000 diluted). Predicted band size: 260KDa. Observed



band size: 260KDa.

Immunoprecipitation analysis of K562 cell lysates using CHD4 mouse mAb.