

APPBP1 antibody



Catalog Number: 108040

Product name

APPBP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

APPBP1 Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue, human brain tissue. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of APPBP1 recombinant protein (Accession Number: NM_003905). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:108040(NAE1 antibody) at dilution of 1:400



IP Result of anti-NAE1 (IP:Catalog No:108040, 4ug; Detection:Catalog No:108040 1:500) with mouse kidney tissue lysate 4000ug.