

FGFR1OP antibody



Catalog Number: 110640

Product name

FGFR1OP antibody

Immunogen

[Human FGFR1OP Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FGFR1OP Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Positive IHC detected in human lymphoma tissue. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of FGFR1OP recombinant protein (Accession Number: NM_194429). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110640(FGFR1OP antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:110640(FGFR1OP antibody) at dilution of 1:50 (under 10x lens)