

transgelin/SM22 antibody



Catalog Number: 107640

Product name

transgelin/SM22 antibody

Recommended Dilution:

WB: 1:500-1:5000

Immunogen

[Human TAGLN Recombinant protein \(GST tag & His tag\)](#)

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

transgelin/SM22 Mouse Monoclonal antibody. Positive IHC detected in human colon tissue, human heart tissue. Positive WB detected in Human colon, human stomach tissue. Observed molecular weight by Western-blot: 18+22 kDa

Preparation

This antibody was obtained by immunization of transgelin/SM22 recombinant protein (Accession Number: NM_003186). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type


Mouse IgG1

Applications


ELISA, WB, IHC

Dilutions


Validations



human colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:107640(TAGLN antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:107640(TAGLN Antibody) at dilution of 1:200 (under 10x lens)



transgelin/SM22 antibody



Catalog Number: 107640

Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog

No:107640(TAGLN Antibody) at dilution of 1:200 (under 40x lens)