

ZGLP1 antibody



Catalog Number: 117050

Product name

ZGLP1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZGLP1 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tissue. Positive WB detected in mouse ovary tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001103167). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



Immunohistochemical of paraffin-embedded human ovary using Catalog No:117050(GLP1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human ovary using Catalog No:117050(GLP1 antibody) at dilution of 1:50 (under 40x lens)



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:117050(ZGLP1 Antibody) at dilution of 1:300