

# FADS3 antibody



Catalog Number: 110450

## Product name

FADS3 antibody

## Immunogen

[Human FADS3 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

FADS3 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, human brain tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 49-51kd

## Preparation

This antibody was obtained by immunization of FADS3 recombinant protein (Accession Number: NM\_021727). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:110450(FADS3 antibody) at dilution of 1:300