

Product name

NFIA antibody

IP: 1:200-1:2000

IHC: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NFIA Rabbit Polyclonal antibody. Positive IP detected in A431 cells. Positive WB detected in A431 cells, HeLa cells, Jurkat cells, L02 cells, mouse liver tissue. Positive IHC detected in human prostate cancer tissue. Observed molecular weight by Western-blot: 60-70kd

Preparation

This antibody was obtained by immunization of NFIA recombinant protein (Accession Number: NM_005595). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:113150(NFIA antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:113150(NFIA antibody) at dilution of 1:25 (under 10x lens)



IP Result of anti-NFIA (IP:Catalog No:113150,

NFIA antibody



Catalog Number: 113150

4ug; Detection:Catalog No:113150 1:500) with A431 cells lysate 2000ug.