

# UFC1 antibody



Catalog Number: 116550

## Product name

UFC1 antibody

Recommended Dilution:

WB: 1:500-1:5000

## Immunogen

IHC: 1:20-1:200

[Human UFC1 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UFC1 Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue, mouse kidney tissue, mouse liver tissue, rat liver tissue. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 21kd

## Preparation

This antibody was obtained by immunization of UFC1 recombinant protein (Accession Number: XM\_017001451). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions


## Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:116550(UFC1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human liver using Catalog No:116550(UFC1 antibody) at dilution of 1:100 (under 10x lens)



# UFC1 antibody



Catalog Number: 116550

---

Immunohistochemical of paraffin-embedded

human liver using Catalog No:116550(UFC1 antibody) at dilution of 1:100 (under 40x lens)