

# ITK antibody



Catalog Number: 111950

## Product name

ITK antibody

## Immunogen

[Human ITK Recombinant protein \(GST tag & His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ITK Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse spleen tissue. Observed molecular weight by Western-blot: 68-72kd

## Preparation

This antibody was obtained by immunization of ITK recombinant protein (Accession Number: NM\_005546). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111950(ITK antibody) at dilution of 1:300