

DNAJC5B antibody



Catalog Number: 109950

Product name

DNAJC5B antibody

IP: 1:500-1:5000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DNAJC5B Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 22kd

Preparation

This antibody was obtained by immunization of DNAJC5B recombinant protein (Accession Number: NM_033105). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109950(DNAJC5B antibody) at dilution of 1:1500



Immunofluorescent analysis of HepG2 cells, using DNAJC5B antibody Catalog No:109950 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



DNAJC5B antibody



Catalog Number: 109950

IP Result of anti-DNAJC5B (IP:Catalog No:109950, 3ug; Detection:Catalog No:109950 1:1500) with mouse testis tissue lysate 4000ug.