

# CD70 Polyclonal Antibody



Catalog Number: 163350

---

**Product name**

CD70 Polyclonal Antibody

**Specificity**

Mouse,Rat

**Antibody description**

Polyclonal antibody to CD70

**Preparation**

Antigen: A synthetic peptide of human CD70.

**Formulation**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

**Clonality**

Polyclonal

**Ig Type**


Rabbit IgG

**Applications**

IF

**Dilutions**

IF 1:50 - 1:200

**Validations**

Immunofluorescence - CD70 Polyclonal Antibody

Immunofluorescence analysis of rat spleen using CD70 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.