

CD164 antibody



Catalog Number: 109060

Product name

CD164 antibody

Specificity

Human; other species not tested.

Antibody description

CD164 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive FC detected in K-562 cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 82 kDa

Preparation

This antibody was obtained by immunization of CD164 recombinant protein (Accession Number: NM_006016). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, WB, FC

Dilutions


Recommended Dilution:

WB: N/AIHC: 1:20-1:200


Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109060(CD164 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:109060(CD164 antibody) at dilution of 1:50 (under 10x lens)



1X10⁶ K-562 cells were stained with 0.2ug CD164 antibody (Catalog No:109060, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

CD164 antibody



Catalog Number: 109060

with dilution 1:1000.