

# SSRP1 Polyclonal Antibody



Catalog Number: 164860

## Product name

SSRP1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to SSRP1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SSRP1 (NP\_003137.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:10 - 1:100

## Validations



Western blot - SSRP1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SSRP1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunohistochemistry - SSRP1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse liver using SSRP1 antibody at dilution of 1:200 (40x lens).

# SSRP1 Polyclonal Antibody



Catalog Number: 164860

---



Immunofluorescence - SSRP1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using SSRP1 antibody .