

Anti-Tnni3 antibody



Catalog Number: 176360

Product name

Anti-Tnni3 antibody

Specificity

Mouse, Rat

Antibody description

Rabbit polyclonal antibody to Tnni3

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from mouse tnni3. aa range:5-54

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1:: Western Blot analysis of mouse cells using Troponin I-C Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000