

FZD5 antibody



Catalog Number: 110760

Product name

FZD5 antibody

Specificity

Human; other species not tested.

Antibody description

FZD5 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, HeLa cells, Jurkat cells, PC-3 cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of FZD5 recombinant protein (Accession Number: BC114346). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:110760(FZD5 Antibody) at dilution of 1:600