

Anti-SLC27A4/FATP4 Rabbit antibody



Catalog Number: 168760

Product name

WB, ICC/IF, IP

Anti-SLC27A4/FATP4 Rabbit antibody

Dilutions

Specificity

WB: 1/1000

Human, Mouse, Rat

ICC/IF: 1/20

Antibody description


IP: 1/20

SLC27A4/FATP4 Rabbit monoclonal antibody

Validations

Preparation

Antigen: Recombinant protein of human SLC27A4/FATP4




Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Western blot detection of SLC27A4/FATP4 in K562, HeLa cell lysates using SLC27A4/FATP4 Rabbit mAb(1:1000 diluted). Predicted band size: 72kDa. Observed band size: 72kDa.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.



Clonality

Monoclonal

Ig Type

Rabbit IgG

Immunocytochemistry of SLC27A4/FATP4(green) in HeLa cells using SLC27A4/FATP4 Rabbit mAb at dilution 1/200, and DAPI(blue)

Applications