

CCNL1 antibody



Catalog Number: 109070

Product name

Polyclonal

CCNL1 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

CCNL1 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, MCF7 cells, PC-3 cells. Observed molecular weight by Western-blot: 46kd

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of CCNL1 recombinant protein (Accession Number: NM_020307). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:109070(CCNL1 antibody) at dilution of 1:400