

Anti-CAP1 Rabbit antibody



Catalog Number: 167470

Product name

Anti-CAP1 Rabbit antibody

Specificity

Human

Antibody description

CAP1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human CAP1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/2000-1/10000

IHC: 1/20-1/100


ICC/IF: 1/50-1/200

IP: 1/20-1/50

Validations



Immunocytochemistry of CAP1(green) in Hela cells using CAP1 Rabbit pAb at dilution 1/50, and



Western blot detection of CAP1 in 293T,MCF7 cell lysates using CAP1 Rabbit pAb(1:1000 diluted).Predicted band size:52KDa.Observed band size:56KDa.