

WDR48 antibody



Catalog Number: 116870

Product name

WDR48 antibody

Immunogen

[Human WDR48 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

WDR48 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in Jurkat cells, HEK-293 cells. Observed molecular weight by Western-blot: 70-80 kDa

Preparation

This antibody was obtained by immunization of WDR48 recombinant protein (Accession Number: NM_001346227). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116870(WDR48 Antibody) at dilution of 1:600



IP Result of anti-WDR48 (IP:Catalog No:116870, 3ug; Detection:Catalog No:116870 1:500) with HEK-293 cells lysate 3000ug.