

GABARAP antibody



Catalog Number: 110770

Product name

GABARAP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GABARAP Rabbit Polyclonal antibody. Positive IF detected in Starvation treated HepG2 cells. Positive WB detected in human placenta tissue, A549 cells, HeLa cells, human heart tissue, mouse brain tissue, mouse liver tissue, mouse pancreas tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 14 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_007278). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110770(GABARAP antibody) at dilution of



1:500

Immunofluorescent analysis of (-20oc Ethanol) fixed Starvation treated HepG2 cells using Catalog No:110770(GABARAP Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)