

# Anti-phospho-EGFR (Tyr1173) Rabbit antibody



Catalog Number: 167770

## Product name

Anti-phospho-EGFR (Tyr1173) Rabbit antibody

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Specificity

Human

## Dilutions

WB: 1/1000

## Antibody description

phospho-EGFR (Tyr1173) Rabbit monoclonal antibody

IHC: 1/20

ICC/IF: 1/20

## Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Tyr1173 of human EGFR

FC: 1/20

IP: 1/20

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage


Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG



Western blot detection of phospho-EGFR (Tyr1173) in HeLa,A549 cell lysates using phospho-EGFR (Tyr1173) Rabbit mAb(1:1000 diluted).Predicted band size:134kDa.Observed band size:175kDa.