

ESRRG antibody



Catalog Number: 110480

Product name

ESRRG antibody

Immunogen

[Human ESRRG Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ESRRG Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue, human kidney tissue. Observed molecular weight by Western-blot: 51kd

Preparation

This antibody was obtained by immunization of ESRRG recombinant protein (Accession Number: XM_017000632). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:110480(ESRRG antibody) at dilution of 1:1000