

# TTL antibody



Catalog Number: 116480

## Product name

TTL antibody

## Immunogen

[Human TTL Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TTL Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, HeLa cells, human brain tissue, mouse brain tissue, rat brain tissue, Y79 cells. Observed molecular weight by Western-blot: 43-45 kDa

## Preparation

This antibody was obtained by immunization of TTL recombinant protein (Accession Number: NM\_153712). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:116480(TTL antibody) at dilution of 1:1000