

TMEM24 antibody



Catalog Number: 116180

Product name

TMEM24 antibody

Immunogen

[Human C2CD2L Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

TMEM24 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 80-100 kDa

Preparation

This antibody was obtained by immunization of TMEM24 recombinant protein (Accession Number: NM_001382613). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116180(C2CD2L Antibody) at dilution of 1:1000