

c-SRC antibody



Catalog Number: 107180

Product name

c-SRC antibody

IHC: 1:20-1:200

IF: 1:50-1:500

Specificity

Human, hamster; other species not tested.

Antibody description

c-SRC Mouse Monoclonal antibody. Positive WB detected in CHO cells. Positive IHC detected in human breast cancer tissue. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of c-SRC recombinant protein (Accession Number: XM_047440405). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



CHO cells were subjected to SDS PAGE followed by western blot with Catalog No:107180(c-SRC Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:107180(c-SRC Antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemistry of paraffin-embedded

c-SRC antibody



Catalog Number: 107180

human breast cancer tissue slide using Catalog No:107180(c-SRC Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (-20oc Ethanol) fixed MCF-7 cells using Catalog No:107180(c-SRC Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)