

ZNF470 antibody



Catalog Number: 117180

Product name

ZNF470 antibody

Immunogen

[Human ZNF470 Recombinant protein \(His tag\)](#)

Specificity

Human, Rat; other species not tested.

Antibody description

ZNF470 Rabbit Polyclonal antibody. Positive WB detected in rat testis tissue. Observed molecular weight by Western-blot: 73-83 kDa

Preparation

This antibody was obtained by immunization of ZNF470 recombinant protein (Accession Number: NM_001001668). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



rat testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:117180(ZNF470 Antibody) at dilution of 1:1000