

CXCR5 Polyclonal Antibody



Catalog Number: 161980

Product name

CXCR5 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to CXCR5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 193-372 of human CXCR5 (NP_001707.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - CXCR5 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using CXCR5 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.