

# XAB2 antibody



Catalog Number: 116980

## Product name

XAB2 antibody

## Immunogen

[Human XAB2 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

XAB2 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HeLa cells. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 100 kDa

## Preparation

This antibody was obtained by immunization of XAB2 recombinant protein (Accession Number: NM\_020196). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116980(XAB2 antibody) at dilution of 1:200



IP Result of anti-XAB2 (IP:Catalog No:116980, 4ug; Detection:Catalog No:116980 1:200) with K-562 cells lysate 3200ug.