

ALG13 antibody

Catalog Number: 107980

Product name

ALG13 antibody

WB: 1:1000-1:10000

IHC: 1:20-1:200

Immunogen

[Human ALG13 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ALG13 Rabbit Polyclonal antibody. Positive IHC detected in human stomach. Positive WB detected in HepG2 cells, Jurkat cells, K-562 cells. Observed molecular weight by Western-blot: 130-135kd

Preparation

This antibody was obtained by immunization of ALG13 recombinant protein (Accession Number: NM_018466). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IF, IHC


Dilutions

Recommended Dilution:


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107980(ALG13 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human stomach using Catalog No:107980(ALG13 antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human stomach using Catalog No:107980(ALG13 antibody) at dilution of 1:100 (under 10x lens)

ALG13 antibody

Catalog Number: 107980

antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using ALG13 antibody Catalog No:107980 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).