

Anti-S6K1 Rabbit antibody



Catalog Number: 168680

Product name

Anti-S6K1 Rabbit antibody

Dilutions

1:1000-1:5000

Specificity

1:50

Human, Mouse, Rat

1:20-1:100

Antibody description

1:20

S6K1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human S6K1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Rabbit IgG


Applications

WB, ICC/IF, FC, IP

Validations



Immunocytochemistry of S6K1(green) in HeLa cells using S6K1 Rabbit mAb at dilution 1/200,



and DAPI(blue)

Western blot detection of S6K1 in Rat Brain lysates using S6K1 Rabbit mAb(1:1000 diluted). Predicted band size: 59kDa. Observed band size: 70kDa.