

PAGE2 antibody



Catalog Number: 113490

Product name

PAGE2 antibody

Immunogen

[Human PAGE2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

PAGE2 Rabbit Polyclonal antibody. Positive IHC detected in human prostate hyperplasia tissue.

Preparation

This antibody was obtained by immunization of PAGE2 recombinant protein (Accession Number: NM_207339). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:
IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human prostate hyperplasia tissue slide using Catalog No:113490(PAGE2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human prostate hyperplasia tissue slide using Catalog No:113490(PAGE2 Antibody) at dilution of 1:50 (under 40x lens)