

# Anti-LAMP2 Rabbit antibody



Catalog Number: 167490

## Product name

Anti-LAMP2 Rabbit antibody

## Specificity

Human

## Antibody description

LAMP2 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human LAMP2

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions


WB: 1/1000

IHC: 1/200

## Validations



Immunohistochemistry of LAMP2 in paraffin-embedded Human lung cancer tissue using



LAMP2 Rabbit pAb at dilution 1/100

Western blot detection of LAMP2 in THP-1 cell lysates using LAMP2 Rabbit pAb(1:1000 diluted). Predicted band size: 45 kDa. Observed band size: 110 kDa.