

BFSP2 antibody



Catalog Number: 108490

Product name

BFSP2 antibody

Immunogen

[Human BFSP2 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

BFSP2 Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue. Observed molecular weight by Western-blot: 45-49 kDa

Preparation

This antibody was obtained by immunization of BFSP2 recombinant protein (Accession Number: NM_003571). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:117135(BFSP2 Antibody) at dilution of 1:3000