

Mouse Vegfd/Figf (His tag) recombinant protein



Catalog Number: 518490

General Information

Gene Name Synonym

Vascular endothelial growth factor D, c-Fos-induced growth factor, FIGF, VEGFD,

Protein Construction

A DNA sequence encoding the Mouse Vegfd/Figf Phe98-Ser206 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Vegfd/Figf (His tag) recombinant protein