

TAF9 antibody



Catalog Number: 115990

Product name

TAF9 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TAF9 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, human colon tissue, human heart tissue, human skeletal muscle tissue, MCF7 cells. Positive IP detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 24 kDa

Preparation

This antibody was obtained by immunization of TAF9 recombinant protein (Accession Number: BC007426). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115990(TAF9 Antibody) at dilution of 1:1000



IP Result of anti-TAF9 (IP:Catalog No:115990, 4ug; Detection:Catalog No:115990 1:1000) with mouse skeletal muscle tissue lysate 2200ug.