

DLG1 antibody



Catalog Number: 109990

Product name

DLG1 antibody

IHC: 1:20-1:200

IF: 1:50-1:500

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DLG1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IHC detected in human heart tissue. Positive WB detected in mouse lung tissue, HeLa cells. Observed molecular weight by Western-blot: 120 kDa or 130-140 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004087).
Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:109990(SAP 97 antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:109990(SAP 97 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (-20oc Ethanol)

DLG1 antibody



Catalog Number: 109990

fixed HeLa cells using Catalog No:109990(SAP 97 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)