

Anti-beta Arrestin 1 (6A1) Mouse antibody



Catalog Number: 168390

Product name

Anti-beta Arrestin 1 (6A1) Mouse antibody

Ig Type

Mouse IgG2b

Specificity

Human

Applications

WB

Antibody description

beta Arrestin 1 (6A1) Mouse monoclonal antibody

Dilutions

WB: 1/500

Preparation

Antigen: Recombinant human β -Arrestin 1 protein.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of β -Arrestin 1 in K562 and Hela cell lysates using β -Arrestin 1 mouse mAb(dilution 1:500). Predicted band size: 51kDa. Observed band size: 51kDa.